


Anatomage[®] TABLET



Explore real anatomy and access terminology instantly. Students can select which structures to display, cut through the anatomy and even open CT or MRI scans for deeper learning.

Content

ANATOMAGE BODIES

1 male, 1 female real life-size cadavers for gross anatomy

PREMADE MATERIALS

74 premade presets and 71 curriculum view sequences

SEGMENTATION

Highly detailed segmentation of 3000+ structures

CURRICULUM

500+ presets of anatomical structures available for supplementing curriculum



CT/MRI SCANS

Applications



BLENDED LEARNING



TEXTBOOK SUPPLEMENT



LECTURE AIDS



HOMEWORK RESOURCES

Instructional & Clinical Tools



VISIBILITY PREFERENCES
FOR SYSTEMS AND
STRUCTURES



ANNOTATION
AND NOTES TOOLS



STRUCTURE
PRONUNCIATION TOOL



DICOM CAPABILITIES



PACS SUPPORTED



SELF-ASSESSMENT AND
QUIZ TOOLS



ACCESS TO ANATOMAGE
SHARE

Hardware

TABLET

PRODUCT DIMENSIONS

| Length: 33 cm (12.85")

| Height: 0.6 cm (0.22")

| Width: 21 cm (8.21")

WEIGHT


| 732g (1.6 pounds)

DISPLAY SIZE

| 14.6"

| I/O

| USB-C



ANATOMAGE EUROPE

Via Chiaravalle 7, 20122 Milan, Italy

+39 02 35 99 05 04

table.europe@anatomage.com

eu.anatomage.com

Anatomage[®]